


Certificate of Analysis

MT-1751

Indoxyl sulfate, potassium salt, [³H]-



Lot #: 394-004-0004-A-20130312-TNG


***Molar Activity:** 0.4 Ci/mmol

Concentration: 1.0 mCi/ml; 628.32 µg/ml

Packaged in: Ethanol : water (1 : 1) solution

Date of Analysis: May 7, 2021

***Radiochemical Purity by HPLC:** 99%



HPLC ANALYSIS LOT 394-004-0004-A-20130312-TNG
File Name: intj731 Date and Time: 5/7/2021 4:24:04 PM
Unit 19 - Radio

Peak #	Area %	Time	Area
1	0.50	4.05000	96.28570
2	0.11	4.50330	20.85948
3	0.24	4.85500	47.13290
4	99.15	5.47670	19079.29495

Totals 100.00 19243.57303

Storage Recommendation: Store at -20°C.

***Quality Level Disclaimer:** Results were obtained utilizing non-qualified instrumentation, methods that have not been validated or transferred, and analyses not performed according to a formal, documented, analytical protocol.

When handling this radiochemical, please refer to the safe handling of isotopes document provided with your shipment. If the document is misplaced, a link to a copy is available at www.moravek.com/customer-service.

In addition to the hazard connected with radiation, this compound also has chemical safety characteristics that need to be taken into account before it can be safely worked with. For assistance with designing a safe working environment and procedures for this and all chemicals at your workplace, please contact OSHA at 800-321-6742(OSHA) or https://www.osha.gov/html/Feed_Back.html.

Caution: This material was manufactured under research grade conditions and is intended for investigational or manufacturing use only. It is pharmaceutically unrefined and is not suitable for use in humans or clinical diagnosis. Responsibility for its use and compliance with federal law rests solely with the purchaser. If the compound and/or any associated information/documentation derived from its use could be used for any potential future filings, research grade material is not recommended. Expiry date not determined. In the absence of stability data, a purity check is recommended before use.